


# Cell Surface Staining Protocol with MultiDots

## Introduction

MultiDots are polymer-protected nanoparticles that encompass inorganic semiconductor nanocrystals. They have increased sensitivity, high fluorescence intensity, resistance to photobleaching, and improved shelf-life stability. A narrow emission spectrum and distinct spectral signature of the MultiDots eliminates the need for complex spectral compensation. For a high-dimensional assay, this allows for the accommodation of more fluorophores within a panel with reduced spread error. MultiDots ease end-users' workflow, reduces the procedure's duration, and improves the staining quality.



## Reagents and Instruments

- Target cell suspension in PBS or buffer of choice
- MultiDot antibody
- Staining Buffer: PBS + 0.5% BSA+ 2 mM EDTA
- 12x75 mm round bottom polystyrene tubes, 50 mL centrifuge tube
- Centrifuge, Vortex, Flow Cytometer

## Direct Labeling


Sample Preparation: Prepare single-cell suspension in staining buffer with desired concentration of cells.

1. Resuspend  $1E+06$  cells in 200  $\mu$ L of staining buffer in a flow tube.
2. Add appropriately titrated MultiDot antibody to the flow tube.

*Note: 10  $\mu$ L of MultiDot antibody for a million cells is recommended*

3. Pipette briefly to mix the conjugate with the cells.

4. Incubate for 25 minutes at room temperature.  
*Note: If cells are very phagocytic, incubate at 4°C*
5. Wash cells by adding 1 mL of staining buffer to the flow tube, followed by centrifugation at 1800 rpm for 5 minutes.
6. Remove the supernatant completely.
7. Resuspend the pellet in 350  $\mu$ L of staining buffer and analyze in flow cytometer.



## Other Applications

- Flow cytometry
- Fluorescent microscopy
- ELISA
- Immunohistochemistry